

# SUCTION - SOS UNIT / ELECTRIC SUCTION

## THERAPY SOS UNIT



4550

Brand New for 2016, the Therapy SOS (Suction Oxygen System) is a neat fully integrated Oxygen and Suction supply device, which is normally powered from an Oxygen Gas Cylinder. It is ideal for all Patient Transfer and Resuscitation Trolleys, or anywhere within a medical environment where piped oxygen or vacuum is unavailable. The unit can be supplied with a variety of fixtures, and carries a 7 Year Function Warranty. The unit is CE Marked, and fully conforms to British Standard. The standard unit is supplied with a 1.5 Metre Hose, however other lengths are available on request.

### Standard 1.5M Hose Assembly

4550	Rail Clamp Mounted Therapy SOS Unit	250.00
4555	Rail Clamp Mounted Therapy SOS Unit for Huntleigh Trolleys	250.00
4556	Rail Clamp Mounted Therapy SOS Unit with Right Angled 1.5m Hose Assembly	250.00
4559	Therapy SOS Unit with Mounting Bracket suitable for Anetic Aid Trolley (QA3 V2 & V3) & Right Angled 1.5m Hose Assembly	390.00
4560	Pole Clamp Mounted Therapy SOS Unit (Horizontal Pole)	300.00
4565	Pole Clamp Mounted Therapy SOS Unit (Vertical Pole)	300.00
4570	V Dovetail Mounted Therapy SOS Unit	250.00

## ELECTRIC SUCTION

A modern range of Electric Suction Units with a lightweight plastic case  
Features include:

- Fully adjustable vacuum without venting to atmosphere.
  - Hydrophobic Filter giving total contamination protection to the motor.
- Available with High Suction only in the form Mains and Mains/Battery Model.



4260

4260	240V Electric High Suction Unit - Single Jar - Castors & Handle	650.00
4261	240V Electric High Suction Unit - Twin Jars - Castors & Handle	690.00
4262	240V Electric High Suction Unit - Without Jar - Castors & Handle	620.00
4263	240V Electric High Suction Unit - Abbots Bracket - Castors & Handle	630.00
4268	240V Electric High Suction Unit - Twin V Bracket - Castors & Handle	630.00
4264	240V Electric High Suction Unit - Single Jar - Freestanding	625.00
4265	240V Electric High Suction Unit - Twin Jars - Freestanding	660.00
4266	240V Electric High Suction Unit - Without Jar - Freestanding	595.00
4267	240V Electric High Suction Unit - Abbots Bracket - Freestanding	615.00
4269	240V Electric High Suction Unit - Twin V Bracket - Freestanding	615.00

## MAINS/BATTERY SUCTION

4280	Electric/Battery High Suction Unit - Single Jar - Castors & Handle	900.00
4281	Electric/Battery High Suction Unit - Twin Jars - Castors & Handle	940.00
4282	Electric/Battery High Suction Unit - Without Jar - Castors & Handle	868.00
4283	Electric/Battery High Suction Unit - Abbots Brkt - Castors & Handle	878.50
4284	Electric/Battery High Suction Unit - Single Jar - Freestanding	884.00

## NEW FOR OLD ELECTRIC SUCTION UNITS

Any old Electric Suction Unit can be used in exchange for new Electric Suction Units

4211	NFO High Electric Suction Unit - Single Jar - Castors & Handle	552.50
4212	NFO High Electric Suction Unit - Twin Jars - Castors & Handle	586.50
4213	NFO High Electric Suction Unit - Without Jar - Castors & Handle	527.00
4214	NFO High Electric Suction Unit - Abbots Bracket - Castors & Handle	535.50
4216	NFO High Electric Suction Unit - Single Jar - Freestanding	531.25
4218	NFO High Electric/Battery Suction Unit - Single Jar - C & H	765.00



4280

SUCTION  
ACCESSORIES  
SEE PAGE 18