



Therapy Equipment Ltd

FULLY CONTROLLABLE
0-15 LITRES PER MINUTE
OXYGEN

VARIOUS RAIL
MOUNTING OPTIONS

DUAL SCALE
(kPa/mmHG)
VACUUM GAUGE



FULLY CONTROLLABLE
0 > 500mmHG ADJUSTABLE
VACUUM

****MAINTENANCE FREE****

THERAPY SOS UNIT

Unit 1, Cranbourne Ave, Cranbourne Ind. Est, Potters Bar, Herts EN6 3JN Tel: 01707 652270

www.therapyequipment.co.uk

FRONT



TOP



REAR




TECHNICAL DATA

- FULLY INTEGRATED OXYGEN AND SUCTION
- IDEAL FOR RESUSCITATION TROLLEY
- VARIETY OF MOUNTING ALTERNATIVES
- UNITS SUPPLIED WITH **7 YEAR FUNCTION WARRANTY**
- NO ANNUAL STRIP DOWN MAINTENANCE REQUIRED
- FULLY CONFORMS TO BRITISH STANDARD (ISO 10079-3 AND ISO 15002)
- CE MARKED
- DIAL SELECT OXYGEN FLOWRATES—OFF/1/2/3/4/5/6/7/8/10/12/15 L/MIN
- ADJUSTABLE SUCTION LEVELS—0->500mmHG
- SUCTION FLOWRATE—> 15 L/MIN
- VENTURI OXYGEN USAGE—APPROX. 40 L/MIN @ FULL VACUUM
- GAUGE ACCURATE TO +/- 2.5%
- STANDARD UNIT SUPPLIED WITH 1.5M HOSE ASSEMBLY—OTHER LENGTHS AVAILABLE
- COMPATIBLE WITH ALL TYPES OF COLLECTION JARS/LINERS
- PIPELINE PROTECTOR (PART NO. 4902) CAN BE FITTED

THERAPY EQUIPMENT LTD

UNIT 1, CRANBOURNE AVENUE
CRANBOURNE INDUSTRIAL ESTATE
POTTERS BAR, HERTS. EN6 3JN
TEL: 01707 652270 FAX: 01707 652622
E-mail: sales@therapiequipment.co.uk
Website: www.therapiequipment.co.uk



Certificate Number MD 595894

Version 1: August 2015