




THE THERAPY EQUIPMENT LTD

F-18 SUCTION PUMP BATTERY

The F-18B is a powerful compact suction machine used for both Acute and Community use. Suitable for Ambulances, rapid response, and minor surgery.



1

Flexible and suitable for a wide range of applications complete with disposable filter and tube sets

2

Supplied with 1L. sterilisable jar with overflow valve

3

Reliable and easy to use

4

Long-life oil-less piston pump motor



Complete with suction regulator and vacuum gauge.



Light weight, quiet, portable and safe



Can be charged on one 12V car-plug

TECHNICAL DATA

Power Supply Wall':	Input 100 V-240v - 50-60hz
Power Supply Wall':	Output 15vdc 4a - 60w
Power Consumption:	Current 0,54a
Fuse:	Ptc 3a 250v
Lead Battery:	12vcc 1,3a
Battery Life:	100 Min.
Battery Charging Time:	90 Min.
Duty Time:	30 Min. On - 30 Min. Off

*If present

Max Aspiration Flow Direct: Measurement (On Aspiration Vent):	24 Lt/Min.
Max Aspiration:	- 0,75 Bar
Ip:	21
Noise DbA:	61,2 - 76
Min. Size Length Width Height:	30x19x16 cm
Max. Size Length Width Height:	30x19x21 cm
Weight:	3,1kg
Vehicle Charging Cable:	12v (Pin 2,1mm / Ø 5,5mm)

AREAS FOR USE


- Minor Surgery
- Post-tracheotomy
- Hospitals and emergency vehicles
- Veterinary and GP Practices

DESCRIPTION	PART NUMBER
F-18 Suction Pump Battery	10501
Pergo 1 Litre Jar	12000
Pergo 1 Litre Jar Liner	12000-01

OUR PRODUCT RANGE



Medical Regulators



Flowmeters



Suction Equipment



Entonox



Portable Suction Device and Accessories



Sundries



Therapy SOS



High/Low Air-Oxygen Mixer

Therapy Equipment Ltd,
Unit 1, Cranbourne Avenue,
Cranbourne Industrial Estate,
Potters Bar, Herts. EN6 3JN

Tel: 01707 652270 Fax: 01707 652622
E-mail: sales@therapyequipment.co.uk

www.therapyequipment.co.uk



Certificate Number MD 555894



Distributed exclusively by
THERAPY EQUIPMENT LTD

