

# REPAIRS RECONDITIONING

It is the policy of Therapy Equipment to only offer a reconditioning service on our own equipment:

- Mediregs
- Flowregs
- Dialregs
- Entonox Mediregs
- Entonox Demand Valve



Whilst we offer a damage repair service (items within 7 years old), we do not offer a reconditioning service on Suction equipment or Flowmeters. We do not offer a repair service on Suction or Flowmeters over 7 years old, or undated.

All other products sent in for repair, including old style Regulators (S60, S400M and Series 90) and other manufacturer's Suction and Flowmeters will, unless otherwise advised, be treated as New for Old. We reserve the right to adjust the price, if the unit is substantially damaged.

## MEDIREG RANGE

Full Workshop Overhaul after 4-year warranty.

All internal seals and High Pressure Valve Seating are replaced, and the unit returned to the Hospital, with a renewed 4 year warranty. Every effort will be made to replace all damaged parts within the price, however should new Gauges or Metering Drums be required, these will be charged separately.

It is recommended that the Regulator be reconditioned a maximum of 2 times, giving the unit a maximum 12 year life. Units returned for reconditioning beyond 12 years old will not be reconditioned, and a New for Old replacement offered.

Medireg

Dialreg

Flowreg

Entonox Medireg

PLEASE NOTE: ALL MEDIREG RANGE EQUIPMENT RETURNED FOR DAMAGE REPAIR WILL AUTOMATICALLY BE RECONDITIONED.



AVAILABLE IN  
THE UK ONLY



## ENTONOX

Full Workshop Overhaul after warranty.

The Valve Seating and Seals are replaced, and the unit returned to the Hospital, with a renewed warranty. The price includes the reconditioning of the Entonox Demand Valve only, and the replacement Hose Assembly will be charged separately.

It is recommended that the Demand Valve be reconditioned a maximum of 3 times, giving the unit a maximum 12 year life. Units returned for reconditioning beyond 12 years old will not be reconditioned, and a New for Old replacement offered.

Demand Valve (Does not include Hose Replacement)

**PLEASE NOTE: THE DEMAND VALVE SHOULD BE FULLY DECONTAMINATED PRIOR TO RETURNING. ANY UNIT FOUND CONTAMINATED WILL BE RETURNED UN-RECONDITIONED.**



## OLD STYLE FLOWMETERS/DIRECT DIAL FLOWMETERS

We are only able to offer damage repair to Therapy Equipment Flowmeters. It is not possible to recondition a Flowmeter after the 7 Year Warranty. Units received over 7 Years old, or where the date of manufacture is unknown, will not be repaired, and a New for Old replacement offered. All Damage Repairs carry a 1 Year Warranty - Function Warranty remains 7 Years from original date of manufacture.

Therapy Equipment Old Style Flowmeter

Therapy Equipment Twin Old Style Flowmeter

Therapy Equipment Low Flow Flowmeter

Therapy Equipment Direct Dial Flowmeter



## DIAMOND RANGE PRODUCTS

Diamond Range Products are covered by a full 7 Year Accidental Damage and Function Warranty. Any Diamond Range Product found damaged should be returned to Therapy Equipment for a Free of Charge Repair/Replacement. All Damage Repairs carry a 1 Year Warranty - Function Warranty remains 7 Years from original date of manufacture.

Therapy Equipment Diamond Range Flowmeter

Therapy Equipment Diamond Range Suction Controller



## SUCTION CONTROLLERS

We are only able to offer damage repair to Therapy Equipment Suction Controllers. It is not possible to recondition a Suction Controller after the 7 Year Warranty. Units received over 7 Years old, or where the date of manufacture is unknown, will not be repaired, and a New for Old replacement offered. All Damage Repairs carry a 1 Year Warranty - Function Warranty remains 7 Years from original date of manufacture.

Therapy Equipment High Suction Controller

Therapy Equipment Low Suction Controller

Therapy Equipment Injector Suction Controller

Therapy Equipment Emergency Injector Suction Unit

## HIGH PRESSURE HOSE ASSEMBLY RECONDITIONING

We are able to offer a reconditioning service on all hose assemblies where the existing fittings are re-used and the hose renewed.

Hose Refurbishment (1m, 2m, 3m and 4m)

AVAILABLE IN  
THE UK ONLY

