



THE THERAPY EQUIPMENT LTD

F-31 SUCTION PUMP AC/DC WITH LITHIUM BATTERY

The New F-31, with a powerful piston motor, is a portable suction machine for Acute and Community use. It is suitable for emergency and minor surgery.



1

Complete with suction regulator and vacuum gauge

2

Flexible and suitable for a wide range of applications complete with disposable filter and tube sets

3

Supplied with 1L. sterilisable jar with overflow valve

4

Long-life oil-less piston pump motor



Reliable, easy to use, quiet and safe



Can be mounted on trolley with 4 castors



Can be charged on one 12V car-plug

TECHNICAL DATA

Power Supply Wall:	240V - 110 50/60Hz	Max Aspiration:	- 0,80 Bar
Power Consumption:	1,5A	IP:	33
Battery Life:	100 Min.	Noise dBA:	61,2
Battery Charging Time:	90 Min.	Min. Size Length Width Height:	45x15x28 cm
Duty Time:	60 Min. On - 60 Min. Off	Max. Size Length Width Height:	45x15x30 cm
Max Aspiration Flow Direct Measurement (On Aspiration Vent):	40 Lt/Min.	Weight:	3,2kg
		Vehicle Charging Cable:	12v (Pin 2,1mm / Ø 5,5mm)

AREAS FOR USE


- Minor Surgery
- Emergency Surgery
- Hospitals and emergency vehicles
- Veterinary and GP Practices

DESCRIPTION	PART NUMBER
F-31 Suction Pump AC/DC with Lithium Battery	10601
F-31 Suction Pump AC/DC with Lithium Battery With Trolley	10701
Pergo 1 Ltr Jar	12000
Pergo 2 Ltr Jar	12500
Pergo 1 Ltr Liner	12000-01
Pergo 2 Ltr Liner	12500-01

OUR PRODUCT RANGE



Medical Regulators



Flowmeters



Suction Equipment



Entonox



Portable Suction Device and Accessories



Sundries



Therapy SOS



High/Low Air-Oxygen Mixer

Therapy Equipment Ltd,
Unit 1, Cranbourne Avenue,
Cranbourne Industrial Estate,
Potters Bar, Herts. EN6 3JN

Tel: 01707 652270 Fax: 01707 652622

E-mail: sales@therapyequipment.co.uk

www.therapyequipment.co.uk



Certificate Number MD 555894



Distributed exclusively by
THERAPY EQUIPMENT LTD

